

*Proven Excellent  
Clinical Performance*

## RADIANT WARMER WITH FIXED BABY CRADLE MODEL | KAY 300

- Micro processor based servo controlled electronic controller
- Skin and manual control
- Two bright digital displays combination of RED and GREEN continues monitoring of baby temp and set temp, set temp easily switchable to  $^{\circ}\text{C}$  to  $^{\circ}\text{F}$
- Baby tray with high standard mattress, four sides clear acrylic walls with front and both sides collapsible, baby tray tilting facility both head up & head down, clear acrylic base for under surface Phototherapy.
- German make ceramic heater with parabolic reflector to maintain uniform temp across the mattress with swiveling source box easy to take X-ray and procedures



 **Kay & Company**

**Head Office:**

25, Netaji Subhash Marg, New Delhi - 110002. INDIA  
Ph. : +91 11 43195600 - 20 Lines  
Email : info@kay-co.com, sales@kaycoindia.com

[www.kaycoindia.com](http://www.kaycoindia.com)

**Works:**

B-316, Okhla Industrial Area, Phase-1  
New Delhi - 110020. INDIA  
Ph. / Fax : +91 11 41613432 / 33  
Email : info@kay-co.com