



*Proven Excellent  
Clinical Performance*

## KAY 300- Series

- Micro processor based servo controlled electronic controller
- Skin and manual control
- Two bright digital displays combination of RED and GREEN continues monitoring of baby temp and set temp, set temp easily switchable to  $^{\circ}\text{C}$  to  $^{\circ}\text{F}$
- Baby tray with high standard mattress, four sides clear acrylic walls with front and both sides collapsible, baby tray tilting facility both head up & head down, clear acrylic base for under surface Phototherapy.
- German make ceramic heater with parabolic reflector to maintain uniform temp across the mattress with swiveling source box easy to take X-ray and procedures



Radiant Warmer Fixed Baby Cradle and Lower Drawers  
Model : KAY 310



Radiant Warmer with Fixed Baby Cradle  
Model : KAY 300



Radiant Warmer Stand with Detachable Baby Cradle  
Model : KAY 320 & KAY 600



Radiant Warmer Stand  
Model : KAY 320