

**Manufacturers and Exporters of Medical  
and Surgical Equipments and Appliance**



**Kay & Company**

Distinct Products.Outstanding Service

► **GENERAL ITEMS**

- Dressing Drum
- Instrument Tray
- Hospital Hollowwares
- Patient Care Utensil
- Instrument Sterilizer
- Operation Theatre Equipments
- Autoclaves



An ISO 9001:2008 & ISO 13485:2003 Certified Company

[www.kaycoindia.com](http://www.kaycoindia.com)



## Dressing Drum



**Jointed**



**Jointed**



**Deep Drawn Jointless**

### KCDD Dressing Drum

SIZES In Inches D X H	SIZES In MM D X H	Economy (Jointed)	Deluxe (Jointed)
5 X 5	145X125	KE 222	KD 444
6 X 5	165X125	KE 223	KD 445
6 X 6	165X150	KE 224	KD 446
7 X 7	190X176	KE 225	KD 447
8 X 6	210X150	KE 226	KD 448
8 X 8	210X200	KE 227	KD 449
9 X 5	235X125	KE 228	KD 450
9 X 9	235X230	KE 229	KD 451
10 X 9	280X230	KE 230	KD 452
11 X 5	285X125	KE 231	KD 453
11 X 9	285X230	KE 232	KD 454
11 X 18	285X460	KE 233	KD 455
12 X 10	302X255	KE 234	KD 456
13 X 6	355X150	KE 235	KD 457
13 X 9	355X230	KE 236	KD 458
14 X 5	360X125	KE 237	KD 459
14 X 9	360X230	KE 238	KD 460
14 X 10	360X255	KE 239	KD 461
15 X 10	390X255	KE 240	KD 462
15 X 12	390X305	KE 241	KD 463
18 X 12	475X305	KE 242	KD 464

- ◆ Product mentioned as "Economy" are made of light gauge.
- ◆ Product mentioned as "Deluxe" are made of higher gauge.
- ◆ Available in 304 grade S.S. with reasonable minimum order.
- ◆ Available Dia form 5" to 18" as per Requirement with reasonable minimum order.
- ◆ Dressing Drum are measured form hinges to hinges.
- ◆ All Sizes are approximate.

### KCDD Dressing Drum

SIZES In Inches D X H	SIZES In MM D X H	Super Deluxe (Seamless)
6 X 6	165X165	KSD 501
8 X 7	205X180	KSD 502
9 X 5	230X125	KSD 503
9 X 9	230X230	KSD 504
11 X 5	280X125	KSD 505
11 X 9	280X230	KSD 506
12 X 10	310X250	KSD 507
14 X 5	370X130	KSD 508
14 X 9.5	370X240	KSD 509
15 X 12	390X305	KSD 510

- ◆ Product mentioned as "Super Deluxe" are made of heavy gauge.
- ◆ Available in 304 grade S.S. with reasonable minimum order.
- ◆ Dressing Drums are measured from hinges to hinges.
- ◆ All Sizes are approximate.



## Instrument Tray



### KCIT Instrument Tray

SIZES In Inches	SIZES In MM	Economy	Deluxe
8 X 6	200X150X48	KIE 510	KID 021
9 X 5	230X125X50	KIE 511	KID 022
9 X 6	230X155X48	KIE 512	KID 023
10 X 8	255X200X48	KIE 513	KID 024
11 X 7	280X180X48	KIE 514	KID 025
12 X 8	305X205X48	KIE 515	KID 026
12 X 10	305X255X48	KIE 516	KID 027
14 X 10	355X255X48	KIE 517	KID 028
15 X 12	380X305X50	KIE 518	KID 029
18 X 12	455X310X50	KIE 010	KID 030

♦ Product Available in with & without cover \*Details mention below



### KCCT Catheter Tray

SIZES In Inches	SIZES In MM	Economy	Deluxe
8 X 3	205X80X36	KCE 381	
8 X 3.25	205X85X45		KCD 482
8.75 X 3.25	222X82X41		KCD483
16 X 4	405X110X48	KCE 384	KCD 484

♦ Product Available in with & without cover \*Details mention below



### KCST Shallow Tray

SIZES In Inches	SIZES In MM	Deluxe
8 X 6	200X150X25	KSTD 201
9 X 6	230X155X25	KSTD 202
10 X 6	255X200X25	KTSD 203
11 X 7	280X180X25	KSTD 204
2 X 8	300X205X25	KSTD 205
12X10	305X255X25	KSTD 206
14X10	355X255X25	KSTD 207
15X12	380X305X25	KSTD 208
18X12	455X305X25	KSTD 209

- ♦ Shallow Tray height can be made from 20mm to 30mm as per requirement.
- ♦ Product mentioned as "Economy" are made of 0.35mm S.S. sheet.
- ♦ Product mentioned as "Deluxe" are made of 0.5mm J4 grade S.S. sheet.
- ♦ Product mentioned as "Super Deluxe" are made of 0.5mm 304 grade S.S.
- ♦ Product can also be made from 0.6mm S.S. sheet to 0.8mm S.S. sheet as per customer requirement with reasonable minimum order.
- ♦ All Sizes are approximate & measured from inner to inner.

\*Details mention below





## Instrument Tray



### **KCDT** Deep Tray

SIZES In Inches	SIZES In MM	Deluxe
9.5X5X2.5	240X130X63	KDTD 921
12X6X2.75	300X155X68	KDTD 922
12X8X2.5	310X195X63	KDTD 923
15X12X2.25	380X305X60	KDTD 924

- ◆ Product Available in with & without cover
- ◆ All sizes are approximate & measured inner to inner.

\*Details mention below



### **KCBT** Baby Tray without Cover

SIZES In Inches	SIZES In MM	Deluxe
18 X 12 X 2.5	455X310X66	KBTD 700

- ◆ Product Available with Handle also.
- ◆ All sizes are approximate & measured inner to inner.



### **KCMT** Mayo Table Tray

SIZES In Inches	SIZES In MM	Deluxe
18 X 12	480X330X22	KCTD 820
23 X 16	610X422X30	KCTD 821

- ◆ All sizes are approximate & measured outer to outer.
- ◆ Product mentioned as "Economy" are made of 0.35mm S.S. sheet.
- ◆ Product mentioned as "Deluxe" are made of 0.5mm J4 grade S.S. sheet.
- ◆ Product mentioned as "Super Deluxe" are made fo 0.5mm 304 grade S.S.
- ◆ Product can also be made from 0.6mm S.S. sheet to 0.8mm S.S. sheet as per customer requirement with reasonable minimum order.

\*Details mention below



### **KCKD** Kidney Dish

SIZES In Inches	SIZES In MM	Economy	Deluxe
6"	160X77X35	KKDE 581	KKDD 681
6" (Specical)	165X88X38	KKDE 582	KKDD 682
8"	205X105X40	KKDE 583	KKDD 683
10"	245X110X40	KKDE 584	KKDD 684
10 X 6 (Basin)	255X155X40	KKDE 585	KKDD 685
12"	305X140X42	KKDE 586	KKDD 686
12" (Special)	305X140X50	KKDE 587	KKDD 687

- ◆ Product Available in with & without cover
- ◆ Product mentioned here as special has extra depth.
- ◆ All Sizes are approximate & measured outer to outer.

\*Details mention below



## Hospital Hollowware

### KCGP Gallipot



KGPD 41 - Gallipot 2^  
KGPD 42 - Gallipot 4^  
KGPD 43 - Gallipot 6^  
KGPD 44 - Gallipot 8^


Capacity  
60ml  
112ml  
175ml  
225ml

D X H Size  
50 x 30 mm  
60 x 40 mm  
75 x 45 mm  
80 x 50 mm

♦ Also available in 304 Grade

♦ All Sizes & Capacity are approximate

### KCFJ Forceps Jar



KFJE 71 - Forceps Jar 2 x 4  
KFJE 72 - Forceps Jar 2 x 5  
KFJE 73 - Forceps Jar 2 x 6  
KFJE 74 - Forceps Jar 2 x 7  
KFJE 75 - Forceps Jar 2 x 7

D X H Size  
50 x 100 mm  
50 x 125 mm  
50 x 150 mm  
50 x 180 mm  
50 x 180 mm

♦ Also available in 304 Grade

♦ All Sizes & Capacity are approximate

### KCI Irrigator (Douche Cane)



KIAE 61 - Irrigator Al.  
KISE 62 - Irrigator S.S.  
KISE 63 - Irrigator S.S.  
KISE 64 - Irrigator S.S.

Capacity  
1000ml  
1000ml  
1500ml  
2000ml

D X H Size  
105 x 135 mm  
125 x 130 mm  
125 x 145 mm

♦ All Sizes & Capacity are approximate



## Hospital Hollowware

### **KCWB** Basin

**KWBE 91 - Basin 10"**  
**KWBE 92 - Basin 12"**  
**KWBE 93 - Basin 14"**  
**KWBE 94 - Basin 16"**



**Capacity**  
**1850ml**  
**2250ml**  
**4200ml**  
**6500ml**

♦ Also available in 304 Grade

♦ Capacity & Size measured approximate

### **KCVB** Bowl

**KVBE 21 - Bowl 4"**  
**KVBE 22 - Bowl 5"**  
**KVBE 23 - Bowl 6"**  
**KVBE 24 - Bowl 7"**  
**KVBE 25 - Bowl 8"**



**Capacity**  
**230ml**  
**330ml**  
**520ml**  
**1010ml**  
**1530ml**



**Capacity**  
**600 ML**



**Capacity**  
**1200 ML**

♦ Capacity & Size measured approximate

### **KCMJ** Measuring Jar



**KMJE 51 - Measuring Jar 500ml**  
**KMJE 52 - Measuring Jar 1000ml**

♦ Measuring Jar is graduated  
Also Available in Plastic Autoclave

### **KCFC** Feeding Cup



**KMCE 53 - Feeding Cup**

Also Available in Plastic Autoclave

### **KCTJ** Thermometer Jar



**KTJD 54 - Thermometer Jar**  
**30 x 100mm**

Also Available in Plastic Autoclave



## Hospital Hollowware

### PLASTIC

#### **KCCP** Cidex Tray




KCPD 992 - Cidex Tray Plastic

Capacity  
10 Ltr.

Size  
845 x 170 x 150 mm

♦ Outer Dimension approximate

#### **KCCS** Cidex Tray



### STAINLESS STEEL

KCSE 882 - Cidex Tray S.S.

KCSE 883 - Cidex Tray S.S.

Size  
705 x 204 x 150 m  
705 x 204 x 150 mm

♦ Also available in customised size & 304 grade Stainless Steel.

♦ All sizes are approximate

#### **KCFVS** Formalin Vapor Sterilizer



SIZES In Inches	Economy	Deluxe
10 X 8 X 6	KFVS 990	KFVS 980
14 X 8 X 6	KFVS 991	KFVS 981
20 X 8 X 6	KFVS 992	KFVS 982
23 X 8 X 8	KFVS 993	KFVS 983
26 X 8 X 8	KFVS 994	KFVS 984
30 X 8 X 8	KFVS 995	KFVS 985

- ♦ Product is made of Acrylic Sheet.
- ♦ Product mentioned as "Economy" is made of 3mm to 3.5mm sheet as per sizes.
- ♦ Product mentioned as "Deluxe" is made of 5mm sheet.
- ♦ Product mentioned as "Super Deluxe" is made of 6mm sheet.
- ♦ Product 10" & 14" are having two trays.
- ♦ Product 20", 23", 26" & 30 are having three trays.
- ♦ Product available in three colors : White, Blue & Green

♦ Measuring Jar is graduated





## Patient Care Utensil

**KCBPF** Bed Pan Female



KBPFE 333

**KCBPM** Bed Pan Male



KBPME 444



KBPFD 334



KBPMD 445



KBPFS 335



KBPMS 446

	Size :
KBPFE 41 - Bed Pan Female	355 X 290 X 70mm
KBPFD 42 - Bed Pan Female	375 x 310 x 70 mm
KBPFS 43 - Bed Pan Female	350 x 285 x 97 mm

	Size :
KBPME 41 - Bed Pan Female	355 X 290 X 70mm
KBPMD 42 - Bed Pan Female	375 x 310 x 70 mm
KBPMS 43 - Bed Pan Female	350 x 285 x 97 mm

♦ Also Available in 304 Grade

♦ All Sizes are Approximate





## Patient Care Utensil


### **KCUME** Urinal Male



**KUME 780**



**KUMD 781**



**KUMS 782**

KUME 780 - Urinal Male

KUME 781 - Urinal Male

KUME 782 - Urinal Male

Capacity

800 ml

1030 ml

1500 ml

### **KCUFE** Urinal Female



**KUFE 680**



**KUFD 681**



**KUFS 682**

KUFE 680 - Urinal Female

KUFE 681 - Urinal Female

KUFE 682 - Urinal Female

Capacity

800 ml

920 ml

1010 ml

♦ Also Available in 304 Grade

♦ Capacity is Approximate

### **KCSM** Spitting Mug



KMSE 910 - Spitting Mug

KMSE 911 - Spitting Mug



Capacity

450 ml

450 ml

D X H Size :


100 x 50 mm

100 x 50 mm

♦ Sizes and Capacity are Approximate



## Instrument Sterilizer



**Non-Electric**

**Jointed**

### **KCISN** Instrument Sterilizer (Non. Elec.)

SIZES In Inches D X H	SIZES In MM D X H	Economy (Jointed)	Deluxe (Jointed)
8 X 4 X 2.5	205X100X60	KISNE 930	KISND 920
10 X 4 X 3	255X100X75	KISNE 931	KISND 921
12 X 6 X 4	305X150X100	KISNE 932	KISND 922
14 X 6 X 4	355X150X100	KISNE 933	KISND 923
16 X 6 X 4	405X150X100	KISNE 934	KISND 924
17 X 8 X 6	430X200X150	KISNE 935	KISND 925
20 X 8 X 6	510X200X150	KISNE 936	KISND 926

- ◆ Non electric Instrument Sterilizer are fuel heated & having inner perforated tray in it.
- ◆ If without perforated tray it can be use as instrument box.



**Electric**

**Jointed**

### **KCISE** Instrument Sterilizer (Elec.)

SIZES In Inches D X H	SIZES In MM D X H	Economy (Jointed)	Deluxe (Jointed)
8 X 4 X 3.5	205X100X90	KISEE 810	KISED 820
8 X 5 X 4.5	205X125X120	KISEE 811	KISED 821
10 X 5 X 4.5	255X125X120	KISEE 812	KISED 822
12 X 6 X 6	305X155X150	KISEE 813	KISED 823
14 X 6 X 6	355X155X150	KISEE 814	KISED 824
16 X 6 X 6	405X155X150	KISEE 815	KISED 825
16 X 8 X 8	405X205X200	KISEE 816	KISED 826
18 X 8 X 8	460X205X200	KISEE 817	KISED 827
20 X 8 X 8	510X205X200	KISEE 818	KISED 828
24 X 8 X 8	610X205X200	KISEE 819	KISED 829
24 X 10 X 8	610X255X200	KISEE 820	KISED 830
24 X 12 X 8	610X305X200	KISEE 821	KISED 831
24 X 12 X 12	610X305X305	KISEE 822	KISED 832

- ◆ Are fitted with bakelite / good quality plastic handle & four legs.
- ◆ Are electrically fitted with 1Kw. - 2Kw. single phase element of 220V/50Hz. as per sizes & three core wire.
- ◆ Are available with thermostat to control temperature of water from size 12" onwards. (Kindly add "T" with regular code.)
- ◆ Are available with thermostat and tray lifting arrangement from size 12" onwards (Kindly add "TL" with regular code)
- ◆ Are available with auto-cut off element for all sizes (kindly add "A" with regular code)
- ◆ Product mentioned as "Economy" are made of light gauge as per size.
- ◆ Product mentioned as "Deluxe" are made of higher gauge as per size.
- ◆ Available in 304 grade S.S. with reasonable minimum order.
- ◆ All Sizes are approximate.
- ◆ Product can be available in various sizes as per customer requirement.



**Electric With  
Cutoff Device**

**Jointed**



## Instrument Sterilizer



### **KCISS** Instrument Sterilizer (Seamless)

SIZES In Inches D X H	SIZES In MM D X H	Economy	Deluxe
12 X 6 X 5	300X155X130	KISSE 780	KISSD 770
16 X 6 X 5	410X155X130	KISSE 781	KISSD 771
17 X 8 X 6	430X205X150	KISSE 782	KISSD 772
20 X 8 X 6	510X205X150	KISSE 783	KISSD 773

- ◆ Product mentioned as "Economy" is having 304 grade body & J4 grade cover & perforated tray.
- ◆ Product mentioned as "Deluxe" is completely made of 304 grade.
- ◆ Product is completely fitted with thermostat & tray lifting arrangement but can be made available without thermostat & tray lifting arrangement.
- ◆ Product also available with plain immersion element or auto-cut off.
- ◆ Product electrically fitted with 1Kw. - 2Kw. single phase element of 220V/50Hz. as per sizes & three core wire.
- ◆ Lifting arrangement can be made as per jointed instrument sterilizer.

## Electric

## Jointless



### **KBUS** Bowl & Utensil Sterilizer


SIZES In Inches D X H	Economy (Jointed)	Deluxe (Jointed)
24 X 24 X 18	KBUS 669	KBUSD 890
24 X 16 X 12	KBUS 670	KBUSD 889

- ◆ Designed for Rapid Sterilization of Instruments, Kidney Trays, Bed Pans, Basins, Bowls, Utensils etc. Its' Body, Dome Shaped Lid and Perforated Tray made of Heavy Gauge Stainless Steel Sheet 304 Grade. Heavy Casted Brass and Chromium Plated Water Inlet Valve with Overflow Arrangement, Water Outlet Valve, Lid & Tray Operated By Strong Hydraulic Pump Foot Operated. A Sturdy Tubular Steel Welded Stand with Ground Leveling Flanges. Immersion Heater Thermostatic Controlled, Front Panel Fitted with Off / On Switch and Pilot Lights. Suitable to operate on 220 – 230 Volts AC Supply.





## Sterile Container System



### **KSBC** Brush Model - Containers

Container for brush 280 x 285 x 85mm, with wall mounting fixtures, Maximum 20 each brushes can be stored in container. To take the brush, please move the handle to right or left direction under the container.

Colors are available in :

Model No.	Color	Color
KC-051		Black
KC-052		Red
KC-053		Blue
KC-054		Green
KC-055		Yellow
KC-056		Grey




### **KSIC** Mini Implant Model - Container

Sizes in MM	
160 x 70 x 60 mm	KSIC - 057
200 x 145 x 60 mm	KSIC - 058



### **KSIC** Implant Model - Container

Sizes in MM	
300 x 149 x 65 mm	KSIC - 059
500 x 169 x 75 mm	KSIC - 060



### **KSWB** Sterilization Wire Basket

Sterilization Wire Baskets, made of Stainless Steel Wire mesh 304 Grade, 100% rust proof and long life, available in following sizes.

Dimensions	Item Code
500 X 250 X 30mm	KSWB - 071
500 X 250 X 50mm	KSWB - 072
500 X 250 X 70mm	KSWB - 073
500 X 250 X 100mm	KSWB - 074
400 X 200 X 30mm	KSWB - 075
400 X 200 X 50mm	KSWB - 076
400 X 200 X 70mm	KSWB - 077
400 X 200 X 100mm	KSWB - 078
250 X 250 X 30mm	KSWB - 079
250 X 250 X 50mm	KSWB - 080
250 X 250 X 70mm	KSWB - 081
250 X 250 X 100mm	KSWB - 082





## Midwifery Kit Box

### KCMBA

- ◆ Deep Drawn Jointless Aluminium Construction
- ◆ Used to keep contents for normal delivery at home in rural areas.
- ◆ Available empty or with contents as per Unicef Guidelines



Aluminium



Plastic



Vinyl

### KCSB Steel Bucket

#### ST. STEEL

#### ST. STEEL

### KCDJ Dressing Jar



	Capacity	Size
KSBE 551 - Bucket	11 Ltr.	252 x 256 mm
KSBE 552 - Bucket	14 Ltr.	285 x 278 mm
KSBE 553 - Bucket	18 Ltr.	320 x 303 mm



KDJD 651 - Dressing Jar
KDJD 652 - Dressing Jar
KDJD 653 - Dressing Jar
KDJD 654 - Dressing Jar



Capacity	Size
1 Ltr.	100 x 125 mm
2 Ltr.	125 x 150 mm
3 Ltr.	150 x 185 mm
5 Ltr.	185 x 215 mm



## Autoclaves



Electric Autoclave



Non-Electric Autoclave

**KCAA**

Autoclave Six Wing Nut  
Aluminium

SIZES In Inches D X H	Capacity in Liters	Economy	Deluxe
12 X 12	20	KAAE 071	KAAD 171
12 X 15	25	KAAE 072	KAAD 172
14 X 15	35	KAAE 073	KAAD 173

- ♦ Aluminium Body Deep Drawn.
- ♦ Ring and Lid Made of Aluminium Casting.
- ♦ Used for Sterilization of Surgical Instruments and Dressings



## Autoclave Pressure Cooker Type



**Electric Autoclave**



**Non-Electric Autoclave**

### Autoclave Cooker Type Aluminium

**KACA**

SIZES In Inches D X H	Capacity in Liters	Deluxe
9 X 9	9	KACA 960
9 X 11	11	KACA 961
9 X 11.5	13	KACA 962
9 X 12	14	KACA 963
9 X 13	15	KACA 964
12 X 10	18	KACA 965
12 X 11	20	KACA 966
12 X 12	23	KACA 967
12 X 14	26	KACA 968
12 X 16	30	KACA 969
12 X 18	34	KACA 970
12 X 22	40	KACA 971
12 X 24	44	KACA 972
14 X 15	35	KACA 973
14 X 17	40	KACA 974
14 X 19	45	KACA 975
14 X 21	50	KACA 976

- ◆ Made of Aluminium but Steel Finish
- ◆ Also available in Electric / Non Electric (Fuel Based)
- ◆ Also available with Semi-Autoclave (Timer Based)
- ◆ Epoxy Coated Finish available in
- ◆ All Sizes are approximate

#### Special Features

- ◆ Used for sterilization of surgical instruments & dressing
- ◆ Seamless Construction (Jointless) from Extra Thick Aluminium Sheet/S.S. Sheet for better Safety & Durability.
- ◆ Good Quality synthetic rubber Gasket for Longer & better Performance and Hygiene.
- ◆ Weight Type Safety Valve & Pressure Regulating Device for better Safety and Functioning.
- ◆ Coloured - coded Pressure Gauge (0-30 lbs/in<sup>2</sup>) to indicate the safe, sterilization and danger zones.
- ◆ Highly strong and durable Vacuum & steam Release Valves.
- ◆ Handles and valve knobs made of very durable, heat resistant and fine quality bakelite.
- ◆ Heating Element made of heavy copper pipe for longer life.
- ◆ Well designed, powder coated Outer Stand to provide better stability & attractive looks.
- ◆ Well packed in individual boxes.
- ◆ Kindly add 'E' for electric & 'NE' for non electric with desired code.



## Vertical Autoclave



### Working pressure & temperature:

- ♦ The working pressure for sterilization process is 15 psi-17psi.
- ♦ The working temperature for sterilization process is 121°C

### **KCVPA** Vertical High Pressure Autoclave

Sizes In Inches D X H	Capacity in Liters	Economy	Deluxe
12 X 20	40	KVPAE 401	KVPAD 501
16 X 24	79	KVPAE 402	KVPAD 502

### Construction:

- ♦ All stainless steel double wall construction, Lid is made of thick die pressed SS 304 sheet in mirror finish. The unit is fitted with 8 wing nuts to hold the lid at high pressure. Automatic Sterilizing cycle at 121°C.

### Autoclave is equipped with:

- ♦ Steam release valve.
- ♦ Double, excess pressure spring safety valve for extra safety.
- ♦ Pressure gauge.
- ♦ Water drain outlet & water inlet.
- ♦ Heating element 3kw. 220v-50Hz

### Sterilizing Cycle & Features:

- ♦ Sterilizing cycle is fully automatic & controlled by Microchip Controller coupled with sensor which senses the temperature of inner chamber.
- ♦ Alternately displays the inner temperature of the chamber and remaining time of the current cycle.
- ♦ Gives alerts after every minute with single beep sound after the sterilization cycle starts.
- ♦ Gives the facility to increase or decrease the sterilizing time as per the requirement of the end user.
- ♦ On completion of cycle the displays shows 'done' with beep sound for 30 sec
- ♦ Depressurizing & unlocking the unit is manual.



## Stainless Steel Vertical Autoclave

### **KCVAS** Vertical Autoclave Wing Nut Stainless Steel

SIZES In Inches D X H	Capacity in Liters	Economy	Deluxe
12 X 22	40	KVAS 600	KVAS 700
14 X 22	55	KVAS 601	KVAS 701
16 X 24	79	KVAS 602	KVAS 702

- ♦ Stainless Steel Body
- ♦ Ring Made of Aluminium Casting
- ♦ Lid Made of Aluminium Casting with S.S. Plate
- ♦ Available in Electric / Non Electric (Fuel Based)
- ♦ Used for Sterilization of Surgical Instruments and Dressings
- ♦ Also Available in Stainless Steel Lid As per Customer Requirement
- ♦ Also available with low water Cutt of Devices

### **Common Features**

- ♦ Available in Electric / Non Electric (Fuel Devices)
- ♦ Kindly add 'E' for electric & 'NE' for non electric with code.
- ♦ There is a Difference of quality in Accessories in Economy & Deluxe model.
- ♦ Electric Autoclave fitted with heating element of 2 Kw single phase 220 V / 50 Hz.
- ♦ Good quality rubber Gasket.
- ♦ Max Pressure 15 PSI
- ♦ Equipped with Safety Pressure Release, Weight, Steam Release / Control Valve
- ♦ Highly strong and durable Vacuum & Steam Release Valves.
- ♦ Weight Type Safety Valve & Pressure Regulating Devices for better Safety and Functioning.
- ♦ Well Designed Epoxy Coated Outer Stand
- ♦ Also Available in Semi Automatic (Timer Based)
- ♦ All Sizes & Capacity are Approximate





## Steam Autoclaves

### **KCSA** KayClave - 9002

The "KAYCLAVE" is a unique vessel that generates steam under pressure and temperature between 120° to 125° Celsius. On reaching the right temperature the instruments are subject to a saturated steam environment for a pre-specified time which results in safe and guaranteed sterilization. The "KAYCLAVE" is sturdily constructed in rustproof, heavy gauge Aluminium for use in rugged condition. The range of "KAYCLAVE" is available in various models and sizes from 4 liters capacity to 40 liters capacity to suit different applications. All the models are equipped with heat resistant handles and other safety features. Needless to mention, stringent Quality Checks are maintained at all stages of production.

#### Features

- ◆ All Aluminum deep drawn Seamless Construction.
- ◆ Stylish, Light, Compact, Portable Easy to use top loading autoclave.
- ◆ Fast 25 minutes Sterilizing cycle at 121° C.
- ◆ Extended life Heating Element and Lid Gasket.
- ◆ Color-coded Pressure Gauge show internal chamber condition during the cycle.
- ◆ Easily operated anywhere, anytime.
- ◆ The KayClave for Dental application are available in two sizes with a special fitment of rotary timer and epoxy finish.



### **KCSA** KayClave - 9001

#### Features

- ◆ All Aluminum deep drawn Seamless Construction.
- ◆ Light compact portable easy to use top loading autoclave.
- ◆ The KayClave offers **fast & silent** sterilization cycle at 121° C.
- ◆ The sterilization Cycle is Microchip controlled with time adjustment facility.
- ◆ The LED Display shows the internal chamber temperature and sterilizing time.
- ◆ It consumes less power & reduce electricity bills.
- ◆ Equipped with double "Over Pressure" protection system.
- ◆ Extended life Lid Gasket and Heating Element.
- ◆ Color-coded Pressure Gauge show internal chamber condition during the cycle.





## **KCFA** Fogger Machine

### Features

- ◆ Generates ultra fine aerosol fog with greater volume & area coverage capacity.
- ◆ An ideal tool for fogging disinfectants, fumigants, biocides, fungicides, sanitizers & pesticides.
- ◆ Fogger Machine enables user to select optimum droplet size & flow control allowing to generate fine mist (1-3micron), medium fine (3-10 micron) suitable for each type of chemical treatment. Aerial fumigation in operation theatres, Laboratories, clinics, disinfection in food & beverage plants, disinfection in animal houses, pest control in hospitality sectors
- ◆ Fogger Machine consume negligible power and virtually require no maintenance. Provides operation theater fumigator to suit the requirements of each and every medical profession. These operation theater fumigator are used in most of the Indian hospitals and dispensaries.

## Ultra Violet Storage Cabinets



**Ultra Violet Storage Cabinet**



**Compact Ultraviolet Storage Cabinet**

- ◆ Guaranteed Inactivation of all Bacteria and Viruses



## Hospital Hollowware

### KCPS Pressure Stoves

Stove, Kerosene, Single Burner  
Pressure Type  
Capacity Approx. :  
1.2 litre (2-1/2 pint).  
Non-Collapsible,  
Complete with spare.



KCPS - 984 Pressure Stove

### KCPS Pressure Stoves

Stove, Kerosene, Four Burner  
Pressure Type, Extra Heavy Duty  
Tank Capacity :  
2.5 litres (4-1/2 pint).  
All brass construction Overall.



KCPS - 985 Pressure Stove

### KCPBP Bed Pan



KBPF - 280 Bed Pan Female



KCBPH - 282 Bed Pan Female Heavy



KCBPM - 283 Bed Pan Male

### KCKD Kidney Dish

Sizes in Inches	Autoclavable	Disposable
6	KCKDA 110	KCKDD 210
8	KCKDA 111	KCKDD 211
10	KCKDA 112	KCKDD 212
12	KCKDA 113	KCKDD 213

#### Autoclavable



#### Disposable



### KCP Urine Pot



KCPU 901 - Urine Pot Two in One

### KCP Urine Pot



KCPU 902 Urine Pot Male

### KCP Urine Pot



KCPUF 903 Urine Pot Female

### KCF Feotoscope



KCFP - 501 Feotoscope Plastic  
KCFP - 501 A Alluminum

### KCFJ Forcep Jar



KCFJ - 502 Forcep Jar

### KCTJ Thermometer Jar



KCFJ - 502 Forcep Jar

Hospital Hollowwares		Operation & Examination Lamp		Nursery Equipments
	Suction Apparatus		Rubber Product	
Hospital Furniture				



## Kay & Company

Manufacturers and Exporters of Medical and Surgical Equipment and Appliances

### Head Office:

25, Netaji Subhash Marg, New Delhi - 110002. INDIA

Ph. : +91 11 43195600 - 20 Lines

Fax : +91 11 23282678 / 23275086

Email : info@kay-co.com, sales@kaycoindia.com

### Works:

B-316, Okhla Industrial Area, Phase - 1, New Delhi - 110020. INDIA

Fax : +91 11 41613432 / 33

Email : info@kay-co.com

Representative:



**Website:** [www.kaycoindia.com](http://www.kaycoindia.com)