

ISO 9001
ISO 13485
CE



Kay & Company

UROLOGY

ANAESTHESIA
& CRITICAL CARE

GENERAL SURGERY

GASTRO-
ENTROLOGY


TRANSFUSION

TAPES & DRESSING

GYNAECOLOGY
& PAEDIATRICS

Medical Disposables

Catalogue 2022



www.kaycoindia.com
www.medicaldisposabledevices.com



ISO 9001
ISO 13485
CE

The right choice is right here.

COMPANY HISTORY AND PROFILE

For over fifty years medical fraternity the world over has trusted Kay & Company as a reliable source for medical instruments, appliances and hospital furniture. Our products are in use in nearly 30 countries of Europe, Africa, Middle East, South East Asia and the Americas. The company caters to healthcare institutions, and clinical laboratories through a global network of direct importers.

Our products are also used by many international "non profit" non governmental organizations providing humanitarian medical services under mandate from global agencies.

Our range caters to all specialties: Cardiovascular/Thoracic, Urological, Gynaecological/Obstetrical, Orthopaedic, Arthroscopic, ENT, Ophthalmic, Plastic Surgery, Dermal, Neurosurgical, Microsurgical and the list goes on.

CERTIFICATIONS

Kay & Company is an ISO 9001:2008 and ISO 13485:2012 certified organization. We are first Indian company to get the prestigious ISO 13485 certification for medical devices. Several of our product categories have the coveted CE certification. This further re-enforces the confidence of importers and end users to indent for our products.

Purchasing executives directly connected with the health care industry are welcome to contact us directly. Our proximity to the Indian dry ports gives us the advantage of competitive pricing and delivery schedules. For direct importers this gives us the benefit of shipping out "less than container load" shipments.



Index

UROLOGY	2-4
ANAESTHESIA & CRITICAL CARE	5-13
GENERAL SURGERY	14-24
GASTRO-ENTROLOGY	25-27
TRANSFUSION	28-33
TAPES & DRESSING	34-36
GYNAECOLOGY & PAEDIATRICS	37-39
ONCOLOGY	40
BLOOD COLLECTION TUBES	41
PATHOLOGY LAB PRODUCTS	42-43

Kay & Company

Kay & Company reserves the right to discontinue any of the products or change any of the specifications provided in this catalogue without prior notice.

High Flow Foley Balloon Catheter

High Flow 2 Way	KS 1001	Foley Balloon Catheter High Flow
High Flow 3 Way	KS 1002	Grand Débit Foley Globo Sonda
Regular Flow 2 Way	KS 1003	Cateter balão High Flow
Regular Flow 3 Way	KS 1004	Cathéter à ballon High Flow

FEATURES

- Used for short/long term urine drainage
- Made from natural latex rubber
- Silicon elastomer coated smooth surface for atraumatic catheterization
- High strength polymer layer in the middle of the catheter ensures wider inner diameter and hence high flow rate
- Minimises encrustation and subsequent catheter blockage and failure
- Smooth eye, ultra thin highly elastic balloon and hard non-return valve for trouble free inflation and deflation
- Sizes (FR) : 8, 10, 12, 14, 16, 18, 20, 24, 26



TWO WAY	
Size in FR	Balloon Capacity
6	3cc
8 - 10	3cc - 5cc
12 - 14	5cc - 15cc
16 - 26	30cc - 50cc

THREE WAY	
Size in FR	Balloon Capacity
16,18, 20, 22, 24 & 26	30cc - 50cc
14	15cc - 30cc



➤ Sterile, Individually packed in peelable pouch.

STERILE

Silicon Foley Balloon Catheter

2 Way	KS 1005	All Silicon Foley Balloon Catheter
3 Way	KS 1006	Todo Es Silicio Foley Globo Sonda
		Todos cateter de silicone Foley Balloon
		Toutes les cathéter Foley Silicon Balloon

FEATURES

- Used for long term urine drainage
- Made of 100% silicon
- Non-toxic, bio-compatible and extra smooth for maximum patient comfort
- Symmetrical balloon, rounded sealed tip and hard non-return valve for trouble free inflation and deflation
- Sizes (FR) : 8, 10, 12, 14, 16, 18, 20, 22, 24, 26



TWO WAY	
Size in FR	Balloon Capacity
6	3cc
8 - 10	3cc - 5cc
12 - 14	5cc - 15cc
16 - 26	30cc - 50cc

THREE WAY	
Size in FR	Balloon Capacity
16,18, 20, 22, 24 & 26	30cc - 50cc



➤ Sterile, Individually packed in peelable pouch pack.

STERILE

Nelaton / Urethral Catheter

KS 1007

Nelaton / Urethral Cateter
Nélaton / Uréthral Catéter
Nelaton / Urethral Catheter
Nélaton / Uréthral Cathéter

FEATURES

- Used for short term bladder catheterization
- Made from non-toxic, non-irritant, medical grade PVC
- Radio opaque line provided through out the catheter for X-ray visualization
- Distal end with lateral eyes and proximal end funnel connector, for maximum drainage and easy connection to the urine bag
- Sizes (FR) : 6, 8, 10, 12, 14, 16, 18, 20

VERSION	EFFECTIVE LENGTH	TOTAL LENGTH
Male	370mm	400mm
Female	170mm	200mm



› Sterile, Individually packed in peelable blister pack.

PHTHALATE FREE

STERILE



Male External Catheter Penile Sheath

KS 1008

Male External Catheter Penile Sheath
Catéter Externo Masculino Con Funda Pene
Cateter masculino externo peniana Bainha
Homme externe cathéter gaine pénienne

FEATURES

- Used for urine drainage in incontinent male patients
- Made from natural latex rubber
- Provided with adhesive strip for proper fixing of the catheter

SIZES	
Small	20mm
Medium	25mm
Large	30mm
Extra Large	35mm

› Sterile, Individually packed in peelable pouch pack.

STERILE



Red Rubber Nelaton Catheter

KS 1010

Red Karet Nelaton Kateter
Catéter de Nelaton de caucho rojo
Cateter vermelho de borracha Nelaton
Cathéter de Nelaton en caoutchouc rouge

FEATURES

- Red natural latex rubber catheter used across the spectrum of health care for all patient types
- High quality construction ensure consistency of use and patient comfort
- Smooth tip and funnel end made from a radiopaque material
- Two large opposing eyes for optimum drainage
- Suitable for male or female catheterization
- Product type - Intermittent Catheter
- Product style - Straight tip
- Sizes (FR) : 8, 10, 12, 14, 16, 18, 19, 20



› Sterile, Individually packed in peelable blister pack.

STERILE

Siliconised Endotracheal Tube - Cuffed & Plain

(Nasal/Oral)


Cuffed | KS 1011 SC

Plain | KS 1012 SP

The Siliconised Endotracheal Tube hidung / lisan
 Tubo Endotraqueal Combinado Con Silicona - Cuffed & Plain
 Endotracheal silikonisada Tube - Cuffed & Plain
 Endotrachéale siliconé Tube - revers & Plaine

FEATURES

- Used to maintain a patent airway and mechanically ventilate patients
- Made from a special non-toxic Silicon and PVC blend
- PHTHALATE FREE - minimizing patient exposure to carcinogenic Pthallate
- Extra Smooth, non-adherent inside/outside surface for ease in intubation, extubation and suctioning
- Kink resistant, thin-walled, thermo-sensitive tube, softens at body temperature
- Smooth polished murphy eye, distal end and tapered balloon end for atraumatic intubation and extubation
- Kink resistant inflation tube to ensure patient safety during cuff inflation and deflation
- Radio opaque line and markings provided to facilitate identification of tube position
- Universal 15 mm connector at proximal end
- E.T.O Sterile and Pyrogen free. For single patient use
- Cuffed Sizes (mm) : 2.5, 3.0, 3.5, 4.0, 4.5, 5.0, 5.5, 6.0, 6.5, 7.0, 7.5, 8.0, 8.5, 9.0, 9.5, 10
- Plain Sizes (mm) : 2.0, 2.5, 3.0, 3.5, 4.0, 4.5, 5.0, 5.5, 6.0, 6.5, 7.0, 7.5, 8.0, 8.5, 9.0, 9.5, 10



› Sterile, Individually packed in peelable pouch

PHTHALATE FREE

STERILE

Endotracheal Tube - Cuffed & Plain

(Nasal/Oral)


Cuffed | KS 1013 C

Plain | KS 1014 P

Endotracheal Tube - Cuffed hidung / lisan
 Tubo Endotraqueal - Cuffed & Plain
 Tubo endotraqueal - Cuffed & Plain
 Tube endotrachéale - revers & Plaine

FEATURES

- Used to maintain a patent airway and mechanically ventilate patients
- Made from non-toxic, non-irritant medical grade clear PVC
- Kink resistant, thin-walled, thermo-sensitive tube, softens at body temperature
- Smooth inside/outside surface for ease in intubation, extubation and suctioning
- Smooth polished murphy eye, distal end and tapered balloon end for atraumatic intubation and extubation
- Kink resistant inflation tube to ensure patient safety during cuff inflation and deflation
- Radio opaque line and markings provided to facilitate identification of tube position
- Universal 15 mm connector at proximal end
- E.T.O Sterile and Pyrogen free. For single patient use
- Cuffed Sizes (mm) : 2.5, 3.0, 3.5, 4.0, 4.5, 5.0, 5.5, 6.0, 6.5, 7.0, 7.5, 8.0, 8.5, 9.0, 9.5, 10
- Plain Sizes (mm) : 2.0, 2.5, 3.0, 3.5, 4.0, 4.5, 5.0, 5.5, 6.0, 6.5, 7.0, 7.5, 8.0, 8.5, 9.0, 9.5, 10



› Sterile, Individually packed in peelable pouch

PHTHALATE FREE

STERILE


Endotracheal Tube - Subglottic Suction Line

KS 1015

Endotracheal tube - subglottic Suction Baris
 Tubo endotraqueal - subglótica Línea de Succión
 Endotraqueal Tube - subglótica da linha de sucção
 Sonde endotrachéale - Ligne d'aspiration sous-glottique

FEATURES

- Used to maintain a patent airway in long term mechanically ventilated patients
- Made from a special non-toxic medical grade PVC
- Special suction port for removal of subglottic secretions using syringe/suction minimize the risk of VAP (Ventilator Associated Pneumonia)
- Extra smooth, non-adherent inside/outside surface for ease in intubation, extubation and suctioning
- Kink resistant, thin-walled, thermo-sensitive tube. Softens at body temperature
- Smooth polished Murphy eye and distal end for atraumatic intubation and extubation
- Radio opaque line and markings provided to facilitate identification of tube position
- Universal 15 mm connector at proximal end
- E.T.O Sterile and Pyrogen free. For single patient use
- Sizes (mm) : 6.0, 6.5, 7.0, 7.5, 8.0, 8.5



Sterile, Individually packed in peelable pouch / blister packed.

PHTHALATE FREE
STERILE

Reinforced Endotracheal Tube - Cuffed & Plain


(Nasal/Oral)

Cuffed | KS 1016 C
Plain | KS 1017 P

Reinforced Endotracheal Tube- Cuffed & Plain
 Tubo Endotraqueal Reforzado Plano
 Endotraqueal reforçado Tube - Cuffed & Plain
 Endotrachéale renforcé Tube - Revers & Plain

FEATURES

- Used to maintain a patent airway and mechanically ventilate patients
- Made from non-toxic, non-irritant medical grade PVC
- Uniform spiral metal-reinforced tube wall to prevent the tube from kinking
- Reinforcement supports tube position and verification during X-ray
- Accurate and safe positioning ensured by graduation marks and tip-to-tip radiopaque line
- High volume / low pressure cuff ensures efficient low pressure cuff seal (For cuffed tubes)
- Universal 15 mm connector at proximal end
- Recommended use of bite block / Guedel airway
- E.T.O sterile and pyrogen free. For single patient use
- Cuffed Sizes (mm) : 3.0, 3.5, 4.0, 4.5, 5.0, 5.5, 6.0, 6.5, 7.0, 7.5, 8.0, 8.5, 9.0, 9.5.
- Plain Sizes (mm) : 3.0, 3.5, 4.0, 4.5, 5.0, 5.5, 6.0, 6.5, 7.0, 7.5, 8.0, 8.5.



Sterile, Individually packed in peelable pouch / blister packed.

PHTHALATE FREE
STERILE

Preformed Endotracheal Tube - Cuffed & Plain

South & North (Nasal/Oral)

Cuffed | KS 1018 KS & KS 1019 C

Plain | KS 1020 S & KS 1021

Preformed Endotracheal Tube - Cuffed & Plain
 Tubo Endotraqueal Preformada - Con Puño/ Plano
 Endotracheal pré-formada Tube - Cuffed & Plain
 Endotrachéale préformes Tube - Revers & Plain

FEATURES

- Used to maintain a patent airway and mechanically ventilate patients
- Made from non-toxic, non-irritant medical grade PVC
- Preformed tubes reduce the risk of kinking and improve access to the operative field
- Tube can be temporarily straightened for ease of suctioning
- Smooth polished Murphy eye, distal end and tapered balloon end for atraumatic intubation and extubation
- Accurate and safe positioning ensured by graduation marks and tip-to-tip radiopaque line
- High volume / low pressure cuff ensures efficient low pressure cuff seal (For cuffed tubes)
- Universal 15 mm connector at proximal end
- E.T.O sterile and pyrogen free. For single patient use

SIZES - ORAL & NASAL											
ID (mm)	3.0	3.5	4.0	4.5	5.0	5.5	6.0	6.5	7.0	7.5	8.0
O.D. (mm)	4.2	4.8	5.5	6.2	6.9	7.6	8.2	8.9	9.6	10.3	10.9

Sterile, Individually peelable pouch packed.

PHTHALATE FREE

STERILE

Laryngeal Mask Airway (PVC & Silicon)

Disposable

KS 1022 P

Laring Masker Airway (PVC & Silikon)
 Máscara Laringea (PVC y silicio)
 Máscara de Laringe Airway (PVC e silicone)
 Masque Laryngé Voies respiratoires (PVC et silicium)

FEATURES

- Supraglottic, non-invasive airway management device
- Made from non toxic, non irritant medical grade Polymer
- Anatomically engineered mask provided with separate inflation line.
- Provided with universal 15 mm connector.
- Marking provided along the length to facilitate identification of tube position
- Pilot balloon to check inflation & deflation of cuff.
- ETO sterile and pyrogen free
- Sizes : 1, 1.5, 2.0, 2.5, 3.0, 4.0, 5.0



Sterile, Individually packed.

PHTHALATE FREE

STERILE

► Silicone Laryngeal Mask Airway (Double Lumen)

KS 1023 S

Silicone Laring Masker Airway (Lumen Ganda)
Máscara De La Laringe De Siliconay (Doble luz)
Vias aéreas da máscara laríngea de silicone (lúmen duplo)
Masque respiratoire laryngé en silicone (Double lumen)

FEATURES

- Supraglottic, non-invasive airway management device
- Made from medical grade silicon
- Anatomically engineered mask provided with separate inflation line
- Provided with universal 15mm connector
- Marking provided along the length to facilitate identification of tube position
- Pilot balloon to check inflation & deflation of cuff
- Double Lumen Airway provides inbuilt lumen for gastric tubes
- ETO sterile and pyrogen free
- Built-in biting block
- Sizes : 3,4 & 5



► Sterile, Individually pouch packed in carton.

PHTHALATE FREE

STERILE

► Tracheostomy Tube - Cuffed & Uncuffed

Cuffed | KS 1024 C

Uncuffed | KS 1025 P

Trakeostomi Tube - diborgol & Uncuffed
La cánula de traqueostomía - Abofeteado y sin manguito
Traqueostomia - Cuffed & uncuffed
Trachéotomie - revers et sans ballonnet

FEATURES

- Used to maintain a patent airway through tracheostomy and mechanically ventilate patients
- Made from non-toxic, non-irritant, medical grade clear PVC
- Kink resistant, thin-walled, thermo-sensitive tube, softens at body temperature
- Lockable obturator with facility to insert guide wires to facilitate intubation using the widely accepted Seldinger guide wire technique
- Anatomically designed tube and obturator with extra smooth and curved distal end for atraumatic intubation and extubation
- Smooth inside/outside surface for ease in intubation, extubation and suctioning
- Kink resistant inflation tube to ensure patient safety during cuff inflation and deflation
- Radio opaque line and markings provided to facilitate identification of tube position
- Universal 15 mm connector at proximal end
- E.T.O Sterile and Pyrogen free. For single patient use
- Sizes (mm) : 3.5, 4.0, 4.5, 5.0, 5.5, 6.0, 6.5, 7.0, 7.5, 8.0, 8.5, 9.0



► Sterile, Individually blister packed.

PHTHALATE FREE

STERILE


Catheter Mount

KS 1026

Kateter Mount
Montaje del catéter
Montagem de cateter
Support de cathéter

FEATURES

- Compatible with all Breathing and Ventilator circuits.
- Double swivel connector and flexible tubing provides mobility and flexibility to the patient end of the circuit.
- The tubes are made of a special flexible co-polymer and have a minimal compliance with low dead space volume.
- The catheter mount is fitted with a double swivel connector
- Each catheter mount has a 22 mm female / 15 mm male end to connect to the breathing system and 15 mm female end to connect to the endotracheal tube or laryngeal mask.



EO Sterile/ Disposable/ Individual Packed.

STERILE

Oxygen Mask

KS 1027

Oxygen Mask
Máscara De Oxígeno
Máscara de oxigênio
Masque d'oxygène

FEATURES

- Used for efficient oxygen administration to patient
- Made from soft, non-toxic, medical grade PVC
- Elastic strap and aluminum nose clip for proper adjustment of mask
- Multi-channel tube ensures regular supply of oxygen even if the tube is accidentally kinked
- Proximal end of tube provided with universal connector for easy connection with oxygen source

SIZES	
Size	Tube Length
Adult	210cm (7ft)
Child	210cm (7ft)



Non- Sterile, individually packed in a polybag.

STERILE


Nasopharyngeal Airways

KS 1028

Nasofaríngeo Airways
Nasopharyngeal Airways
nasofaríngeas Airways
nasopharyngés Airways

FEATURES

- Used to secure and open airways for patients in case of emergencies
- Made from soft, non-irritant medical grade PVC
- Anatomically designed with smooth round edges, thin wall and kink resistance for maximum patient comfort
- E.T.O. sterile & pyrogen free. For single use only
- Sizes (mm) : 3.0, 3.5, 4.0, 4.5, 5.0, 5.5, 6.0, 6.5, 7.0, 7.5, 8.0, 8.5



› Sterile individually peelable pouch packed.

PHTHALATE FREE

STERILE

Guedel Airways / Oropharyngeal Airways

KS 1029

Guedel Airways / Orofaríngeo Airways
Cánulas De Guedel / Orofaríngeos Airways
Guedel Airways / orofaríngeos Airways
Guedel Airways / oropharyngés Airways

FEATURES

- Used to maintain an unobstructed oro-pharyngeal airway
- Made from rigid, non-toxic medical grade PVC
- Colour coded in-built bite blocks to avoid airway occlusion and biting of tongue

SIZES	Colour Code
000	-
00	Blue
0	Black
1	White
2	Green
3	Yellow
4	Red
5	Orange



› Sterile, individually packed in peelable pouch pack.

PHTHALATE FREE

STERILE

Nasal Oxygen Cannula

Twin Bore Oxygen Delivery Set

KS 1030 A

Nasal Oxygen Cannula Twin Bore Oxygen Delivery Set
Canula Nasal De Oxígeno Entrega Con Doble Juego
Nasal cânula de oxigênio Bore Twin Set entrega de oxigênio
Canule nasale d'oxygène Bore Twin Set oxygène de livraison

FEATURES

- Used for efficient oxygen administration to patients
- Made from soft, non-toxic, medical grade PVC
- Twin nasal tips ensure equal volume of oxygen to both nasal passages for greater patient comfort
- Multi-channel tube ensures regular supply of oxygen even if the tube is accidentally kinked
- Proximal end of tube provided with universal connector for easy connection with oxygen source

SIZES

Size	Tube Length
Adult	210cm (7ft)
Child	210cm (7ft)
Infant	210cm (7ft)



Non-Sterile, individually packed in a polybag.

STERILE

Nebulizer Mask

Child/Adult

KS 1031 A

Nebulizer Masker
Máscara Nebulizadora
Máscara Nebulizador
Masque De Nébulisation

FEATURES

- Specially designed nebulizer chamber with nebulization rate of approximately 3 cc in 10 minutes.
- Average particle size generation is approximately 2 microns.
- Complete kit consists of: Nebulizer mask, Nebulizer chamber and 210 mm Star Lumen Tube.
- Sizes : Adult and paediatric



Sterile, individually packed in blister pack.

STERILE

Suction Catheters


Plain / Normal | KS 1032 P
Thumb Control | KS 1033 A TC
Finger Tip Control | KS 1034 B FT
Atraumatic | KS 1035 T

Suction Catheters
Catéteres De Succión
Cateteres de sucção
Los catéteres de succión

FEATURES

- Used for removal of secretions from the patient's mouth, trachea and bronchial tube
- Made from non-toxic, non-irritant medical grade PVC
- Lateral eye at distal end for effective suctioning
- Proximal end of tube provided with universal connector for easy connection with suction source
- Optional: Thumb controlled, Finger tip controlled, Atraumatic tip
- Available with graduation

VERSION	EFFECTIVE LENGTH	TOTAL LENGTH
Plain	500mm	530mm
TC	480mm	530mm
FT	470mm	530mm
Autromatic	500mm	530mm



› Sterile, individually packed in blister pack.

STERILE

Intubating Stylet

KS 1036
intubasi stylet
intubación estilete
intubating Stylet
stylet d'intubation

FEATURES

- Used along with Endotracheal Tube to aid during intubation
- Smooth plastic coated with flexible malleable aluminum design, allows the stylet to be easily adapted to any desired shape
- Soft distal tip to reduce the possibility of intubation damage
- Smooth surface for ease in introduction and removal from ET tube
- Clean packed

SIZES
6 Fr/Ch (2.0 mm)
10 Fr/Ch (3.3 mm)
14 Fr/Ch (4.7 mm)



› Individually pouch packed.

STERILE

Visionary Bougie Vented

KS 1037

Visioner Bougie Vented
Visionario Bougie ventilado
Visionário Bougie Vented
Visionnaire Bougie ventilé

FEATURES

- Provides rigidity to endotracheal tubes to facilitate endotracheal intubation
- Ultra smooth surface for low friction
- Atraumatic angled tip for ease of intubation and to minimise risk of patient trauma
- Marking provided along the length to facilitate identification of tube position
- Angled tip for ease of intubation
- Vented lumen for oxygen feed
- Sizes (FR) : 6, 10, 15.



› Sterile, individually packaged in easy to open sachets.

STERILE

Thoracic Trocar Catheter

KS 1039

Thoracic Trocar Kateter
Catéter Trocar Torácica
Trocar torácica Cateter
Thoracique trocars catheter

FEATURES

- Used for rapid non-operative plueral and chest drainage
- Made from non-toxic, medical grade PVC
- Catheter is snugly fixed over the trocar, and provided with large smooth eyes for efficient drainage
- Smooth taper shape tip
- Markings at 5, 10, 15 and 20cm from the distal end to measure depth of placement
- Proximal end fitted with taper shape connector for easy connection to drainage bottle/bag

SIZES										
Size	8	10	12	14	16	20	24	28	32	36
Length (mm)	230	230	230	230	230	390	390	390	390	390



› Sterile, individually packed in peelable pouch pack.

PHTHALATE FREE

STERILE

Thoracic / Chest Drainage Catheter


KS 1040

Thoracic / Chest Drainage Catheter
Catéter De Drenaje Torácico
Thoracic / Chest drenagem por cateter
Thoracique / Chest cathéter de drainage

FEATURES

- Used for post-operative drainage after cardio thoracic surgery
- Made from non-toxic, medical grade PVC
- Smooth and round open distal end with large eyes for maximum drainage
- Catheter is marked at every 2 cm from the last eye to ascertain the depth of placement
- Proximal end fitted with taper tongue connector for better forceps grip and smooth penetration
- Optional: Also available in right angle

SIZES								
Size	12	14	16	20	24	28	32	36
Length (mm)	450	450	450	450	450	450	450	450



› Sterile, Individually packed in peelable pouch.

PHTHALATE FREE

STERILE

Peak Flow Meter

With coloured Peak Flow Zones

KS 1042

Pengukur Aliran Puncak
Medidor de flujo máximo
Medidor de fluxo máximo
Débitmètre de pointe

FEATURES

- Peak Flow Meter is easy to use at the doctor's clinic, home or anywhere.
- The Peak flow readings provide a simple, objective parameter to monitor the patient's asthma during visit to the Doctor
- Range: 60-900 Litres / Min (EU Scale)



Breathing Exerciser

KS 1043

Berolahraga Pernapasan
Ejercitador Respiratorio
respiração Exerciser
respiration exerciseur

FEATURES

- Used as a breathing exerciser to increase lung capacity in patients with decreased lung function
- Made from rigid, non-toxic medical grade PVC
- Specially designed with 3 stage chambers to provide more efficient exercise to patient at a step up basis
- Device comprises a base and a central part with a 3 stage chamber, each containing spheres of different sizes, weight and colour
- Provides lung exercise on a step-up basis, as the spheres rise proportionately to the breathing effort exerted by the patient
- Can be dis-assembled into parts for cleaning and disinfection



➤ Individually packed in a polybag.

Urine Collection Bag with Bottom Outlet

KS 1044

Tas Pengumpul Urine
Bolsa de recogida de orina
Saco de coleta de urina
Sac de collecte d'urine

FEATURES

- Soft and kink resistant PVC tubing
- With non-return valve
- Bag graduated in ml to indicate the quantity of urine collected
- Conical inlet connector with cap
- Tube length : 90 & 100 cm
- Available sterile or non sterile.
- For short and long term urine collection
- Reinforced hole for hanging
- Capacity : 2000 ml



STERILE

Urine Collection Bag - Paediatric

KS 1044 P

Tas Pengumpul Urine
Bolsa de recogida de orina
Saco de coleta de urina
Sac de collecte d'urine

FEATURES

- For short term urine collection
- Suitable for both male and female infants
- Capacity : 100 ml
- Provided with special adhesive for better grip and easy removal from soft skin of the infants



STERILE


Urine Collection Bag

KS 1045

Tas Pengumpul Urine
Bolsa de recogida de orina
Saco de coleta de urina
Sac de collecte d'urine

FEATURES

- Soft and kink resistant PVC tubing
- With non-return valve
- Bag graduated in ml to indicate the quantity of urine collected
- Conical inlet connector with cap
- Tube length : 90 & 100 cm
- Available sterile or non sterile.
- For short and long term urine collection
- Reinforced hole for hanging
- Capacity : 2000 ml



Chest Drainage System /Under Water Sealed Drainage System

KS 1046

Dada Sistem Drainase / Under Water Sealed Sistem Drainase
Pecho Sistema de Drenaje / Under Water Sistema de Drenaje Sellado
Peito Sistema de Drenagem / Under Water sistema de drenagem estanque
Système de drainage thoracique / Système de drainage scellés sous l'eau

FEATURES

- Used to remove air (Pneumothorax) or fluid (Pleural Effusion, blood) or pus (empyema) from the intrathoracic space
- Made from non-toxic medical grade PVC
- Bottle with graduation markings to determine precise drain volume
- Clear initial level markings for easy and correct usage
- Suction port to connect bottle with suction
- Universal connectors and kink-resistant tubings for smooth unrestricted flow
- Hangers available for ease of use



SIZES

Single Chamber	500 ml
Single Chamber	1000 ml
Double Chamber	2000 ml

Individually packed in Inner Carton.

Closed Wound Suction Unit

KS 1047

Tertutup Luka Suction Satuan
Unidad De Aspiración Herida Cerrada
Fechado Unidade de sucção da ferida
Fermé Unité d'aspiration des plaies

FEATURES

- Used for closed wound drainage under negative pressure post-surgical patients
- Made from non-toxic, non-irritant PVC
- Radio opaque drainage catheters have multi-perforation for efficient drainage
- Option to clinicians to use one or two catheters simultaneously
- Set contains – Easy to use flexible bottle chamber with hanger, connecting tube with matching size connectors, needle with perforated catheter, one extra separate perforated catheter

SIZES

500ml

600ml

10

Black

12

White

14

Green

16

Orange

18

Red

14

Green



› Sterile, individually packed in peelable pouch pack.

PHTHALATE FREE

STERILE

Yankaur Suction Set

With Standard Handle

KS 1048

Yankaur Suction Set With Standard Handle
Juego De Aspiración Yankaur Con Asa Standard
Yankaur Conjunto de Sucção Com Punho Padrão
Yankaur aspiration avec poignée standard

FEATURES

- Used for removal of secretion and blood pass during surgical procedures
- Made from non-toxic, non-irritant PVC
- Flexible kink resistant tube provides uninterrupted flow during suction
- Soft connectors at both ends of the tube for easy connection
- Available in Standard / Round Tip Handle
- Suitable for use with or without vent control facility

SIZES

Tube Length	2500mm (2.5mtr.)
Tube Diameter	9mm



› Sterile, individually packed in peelable pouch pack.

STERILE

H.I.V Kit / EBOLA Kit

KS 1049

H.I.V Kit / Ebola Kit
H.I.V Kit / Kit EBOLA
Kit H.I.V / EBOLA Kit
Kit H.I.V Kit / EBOLA

FEATURES

- Complete solution for clinicians while examining HIV/Ebola ptients
- Kit consist of
 - Proto Gown -1Pc.
 - Hood Cap -1Pc.
 - Face Mask (with Bacterial Filter) -1Pc.
 - Shoe Legging (Knee Length) -1Pair
 - Latex Gloves -1Pair.
 - PE Sheet -1Pc.
 - Goggles -1Pc.
 - Carry Bag -1Pc.
- Optional: Other customized kits available for Gynecology & Obstetrics, General Surgery, Orthopaedics, Urology



➤ Sterile, Individually packed in peelable pouch.

STERILE

Surgical Gown

KS 1050

Gaun Bedah
Bata quirúrgica
Vestido cirúrgico
Blouse chirurgicale

FEATURES

- Used as a protective barrier to prevent cross-contamination
- Made from soft non-woven fabric with stitched interlocking
- Light weight, with soft elastic cuff and tie straps for better fitting and complete body wrap



KaySafe Surgical Gloves

KS 1051

Sarung tangan bedah
Guantes Quirurgicos
luvas para cirurgia
Gants chirurgicaux

FEATURES

KaySafe Surgical Gloves are designed to be perfectly anatomical for comfortable wrinkle-free fitting. The extremely soft rubber reduces finger fatigue and increases efficiency. The special rubber film used in every pair of KaySafe Gloves adds sensitivity while giving complete protection to the hand. Apart from ergonomic features, our stringent quality control measures make it the safest protective hand wear for surgeons. Lower modulus for better comfort

- Low powder content that minimizes skin irritation
- Thicker gloves for enhanced safety
- Material: Natural Rubber Latex
- Colour: Creamy White
- Design: Hand Specific, Beaded Cuff, Micro rough textured
- Extractable Protein Level: <200 micro gram/ dm².
- Sterilization: Ethylene Oxide / Gamma irradiated as per customer requirement
- Labelling: Shall comply with the corresponding labelling specification and customer requirement
- Shelf Life: 5 years from the date of manufacture.
- Storage condition: Shall be stored in cool dry place and away from direct light



➤ Individually packed

Examination Gloves

KS 1051 E

Sarung tangan pemeriksaan
Luvas de Exame
Guantes de examen
Gants d'examen

FEATURES

- Examination gloves are available in 3 materials Latex, Nitrile (Latex Free) & Vinyl (Latex Free) type. The examination gloves are tested and proved safe against dermal sensitization and skin irritation.
- Material: Natural Rubber Latex
- Colour: Pale Yellow, the color may vary due to storage time and conditions
- Design: Straight Fingers, Thumb and Fingers in one plane, Fits either hand (Ambidextrous) Rolled Rim.
- Powder Content: <10.0 mg/dm² max
- Extractable Protein Level: <200 µg/ dm²
- Labeling: Shall comply with the corresponding labeling specification and customer requirement.
- Shelf Life: 5 years from the date of manufacture.
- Storage condition: Shall be stored in cool dry place and away from direct light

Nitrile Gloves



Latex Gloves

➤ Individually packed

➤ Rubber Gloves / Household Gloves

KS 1053

Sarung tangan karet
 Guantes de goma
 Luvas de borracha
 Des gants en caoutchouc

FEATURES

- Rubber Gloves are used for cleaning work in Hospital and Laboratories,
- It provides excellent abrasion, cut and tear resistance with strong grip and dexterity.
- Wearing gloves protects the hands from harsh detergents and other cleaning products.
- Different length and thickness of the gloves provide various protection.
- The higher thickness & longer cuff gloves provide better protection and are more useful for all chores.
- Available Sizes: Small, Medium & Large



➤ Individually packed

➤ Surgeon Cap / Nurses Cap

Surgeon Cap | KS 1054

Nurses Cap | KS 1055

Gorro De Enfermeras
 Gorro De Cirujano
 Cap cirurgião / tampão das enfermeiras
 Chirurgien CAP / infirmières Cap

FEATURES

- Used to protect against contamination from the head/hair
- Made from non-woven fabric
- Provided with stitched interlocking and elastic strip on the back for proper grip
- Optional: Available in Blue /Green colour



➤ Individually packed.

Surgeon Face Mask

KS 1056

Dokter bedah Masker Wajah
Mascara Facial De Cirujano
Cirurgião Máscara Facial
Visage Chirurgical Masque

FEATURES

- Used as a protective barrier to prevent cross-contamination
- Made from non-woven fabric and available in two-layers and three - layers configuration
- Ultrasonically welded layers for efficient bacterial filtrations
- Provided with a nose clip to prevent fogging of spectacles and to maximize filtration efficiency



➤ Individually packed.

Shoe Cover

KS 1057

Shoe Cover
Cubierta De Zapatos
Shoe Cover
Couvre-chaussure

FEATURES

- Used to protect against contamination from footwear
- Made from non-woven fabric/polyethylene, on a fully automatic ultrasonic sealing machine
- Elastic strip for proper grip



➤ Individually packed.

Surgical Blade

KS 1041

Pisau bedah
Hoja Quirúrgica
Lâmina cirúrgica
Lame chirurgicale



FEATURES

- Disposable Surgical Blades
- Available Size Numbers 10, 10A, 11, 11P, 12, 12B/ 12D, 14, 15, 15T, 15C & 16 fit handle numbers 3,3L, 5,7 & 9.
- Numbers 18, 19, 20, 21, 22, 23, 24 & 25 fit handle numbers 4, 4L & 6 Individually V.C.I paper and foil packed.
- Supplied sterile and non sterile, in box packing of 100

Surgical Blade With Handle

KS 1050

Pisau bedah Dengan Handle
Hoja Quirúrgica Con Asa
Lâmina cirúrgica Com Handle
Lame chirurgicale Avec Poignée

FEATURES

- Disposable, Sterile, ready to use
- Fitted on handle No. 3 or 4



➤ Pack of 20. Carton of 200pcs.

Abdominal Drainage Kit

KS 1052

Perut Drainase Kit
Kit Drenaje Abdominal
Kit de Drenagem abdominal
Kit de drainage abdominale

FEATURES

- Used for abdominal drainage
- Made from non-toxic, medical grade PVC
- Smooth and round distal end with large atraumatic eyes for efficient drainage
- Radio-opaque line provided on the catheter for x-ray visualization
- Marked at every 2cm from the cast eye to ascertain the depth of placement
- Proximal end provided with universal taper connector for easy connection to drainage bag
- 2000ml capacity bag with graduated scale



Box of 10; Master Carton of 100.

Chest Bag / Water Sealed Drainage System

KS 1054

Chest Bag / Water Sealed Drainage System
Bolsa De Pecho/ Sistema De Sellado Pra El Agua Del Drenaje
Saco Peito / Água Selado Sistema de Drenagem
Poitrine Sac / Eau Sealed système de drainage

FEATURES

- Used for collection of wound and empyema
- Made from non-toxic, non-irritant PVC
- 1000ml capacity bag with graduated scale



Individually packed in a polybag.

Nasogastric Tube / Ryle's Tube

Silicon Elastomer Coated

KS 1058 E

Uncoated | KS 1059

Ryle Tabung Silicon Elastomer Coated

Ryle's Tube / Tubo De Ryle Revestimiento de Elastómero de Silicio


Tubo de Ryle Revestimento de Elastômero de Silício

Tube de Ryle Revêtement en élastomère de silicium

FEATURES

- Used for nutritional feeding purposes or aspiration of intestinal secretions
- Made from non-toxic, medical grade PVC
- Enhanced patient safety
- Coned distal end with steel balls sealed into the tube to facilitate easy insertion
- Four lateral eyes for efficient drainage
- Proximal end universal funnel connector
- Markings at 50, 60 and 70 cm and radio opaque line for accurate placement
- Optional: KS 1058 L available with Luer mount fitting
- Ryle's Tube Uncoated KS 1059 also available

EFFECTIVE LENGTH	TOTAL LENGTH
119cm	122cm



› Sterile, Individually packed in peelable pouch & blister pack.

PHTHALATE FREE

STERILE

Salem Sump Tube / Irrigation Tube

Gastro Duodenal Dual Flow Tube

KS 1060

Salem Sump Tube/Irrigation Tube Gastro Duodenal Dual Flow Tube

Salem Sumidero Tubo / Tubo De Irrigación Tubo Gastroduodenal Dual Flow

Salem Sump tubo / tubo de irrigação Gastro Duodenal dupla Fluxo Tubo

Salem puisard Tube / Irrigation Tube Gastro duodénal Dual Flow Tube

FEATURES

- Used to prevent the accumulation of gastric secretions and air
- Made from non-toxic medical grade PVC
- Double lumen tube throughout the length provides airvent when the stomach has been evacuated of accumulated liquid, air or small particulate matter
- Coned distal end with ten lateral eyes for faster drainage
- Universal funnel connector along with double taper connector at proximal end for easy extension
- Markings at 45, 55, 65 and 75 cm and radio opaque line for accurate placement

EFFECTIVE LENGTH	TOTAL LENGTH
119cm	122cm



› Sterile, Individually packed in peelable pouch pack.

PHTHALATE FREE

STERILE

Levin's Tube / Stomach Tube

KS 1061	Levin's Tube / Stomach Tube
	Tubo Levin's / Tubo De Estomago
	Tubo / tubo no estômago de Levin
	Tube / Ventre Tube de Levin

FEATURES

- Used for gastro-intestinal drainage
- Manufactured from non toxic, medical grade PVC
- Coned distal end is coned with four lateral eyes for efficient drainage
- Proximal end universal funnel conector
- Markings at 45, 55, 65 and 75 cm and radio opaque line for accurate placement

EFFECTIVE LENGTH	TOTAL LENGTH
119cm	122cm



6 Green	8 Blue	10 Black	12 White	14 Green	16 Orange	18 Red	20 Yellow
------------	-----------	-------------	-------------	-------------	--------------	-----------	--------------

› Sterile, Individually packed in peelable pouch pack.

PHTHALATE FREE

STERILE


Kehr's "T" Tube

KS 1062	Kehr's "T" Tube
	Tubo De Kehr "T"
	Tube "T" de Kehr
	"T" Tube de Kehr

FEATURES

- Used for short term post operative drainage of common bile duct
- Made from soft non-toxic, medical grade PVC
- Radio opaque line provided throughout the tube for X-ray visualization
- Super smooth tube surface to prevent clogging

SIZES	
Size	10FG , 12FG , 14FG , 16FG & 18FG
Length	20 x 60cm



› Sterile, Individually packed in peelable pouch pack.

PHTHALATE FREE

STERILE

Rectal Catheter

KS 1063

Kateter dubur
Sonda Rectal
cateter rectall
cathéter rectal

FEATURES

- Used for introduction of enema solution into rectum
- Made from non-toxic, non-irritant medical grade PVC
- Open distal end is open with one lateral eye for efficient drainage
- Proximal end with colour coded funnel connector for easy connection to the suction equipment
- Funnel connectors are colour coded for instant size identification
- Optional: Available in 40cm length

EFFECTIVE LENGTH	TOTAL LENGTH
27cm	40cm



› Sterile, Individually packed in peelable pack.

PHTHALATE FREE

STERILE

Infant Feeding Tube Plain / Elastomer Coated

With Graduation | KS 1064 G

Without Graduation | KS 1065


Elastomer Coated | KS 1066 E

Infant Feeding Tube
Tubo De Alimentación Infantil
Infant tubo de alimentação
Infant Feeding Tube

FEATURES

- Used for neonates and paediatric feeding
- Made from non-toxic, medical grade PVC
- Smooth coned distal end for atraumatic intubation, with two lateral eyes for efficient drainage
- Markings at 20cm and radio opaque line for accurate placement and X-ray visualization
- Proximal end with luer mount for connection with feeding device
- Marking from the tip 5 cm to 40 cm

EFFECTIVE LENGTH	TOTAL LENGTH
50cm	52cm



› Sterile, Individually packed in peelable blister pack.

PHTHALATE FREE

STERILE

I.V Cannula With Injection Port

KS 1067

I.V Cannula With Injection Port
Cánula I.V. Con Puerto De Inyección
I.v Cannula com injeção Porto
I.v canule à injection Port

FEATURES

- Used for intra venous infusion of medicinal fluids
- Made from non-toxic medical grade PTFE
- Formulated cannula tip and smooth surface for atraumatic catheterization
- Super sharp needle with two cutting edges
- Provided with injection port with non-return valve for extra medication



SIZES

Size in FG	14	16	18	20	22	24
Colour Code	Orange	Grey	Green	Pink	Blue	Yellow
Length	45mm	45mm	45mm	32mm	25mm	19mm

› Sterile, individually packed in peelable blister pack.

STERILE

Pressure Monitoring Line

KS 1068

Pressure Monitoring Line
Línea De Monitoreo De La Presión
Linha de Monitorização da Pressão
Surveillance de la pression de ligne

FEATURES

- Used for high pressure monitoring, and connection between syringe infusion pump and patient
- Made from non-toxic medical grade PVC
- Can withstand a pressure of 580 PSI
- Provided with universal male luer lock and female luer lock at either end to ensure securely fit all standard equipment
- Optional lengths available



SIZES

Length	110cm	150cm	200cm
Prime Volume	0.9ml	1.40ml	1.80ml

› Sterile packed in peelable pouch pack.

STERILE


➤ Elastic Adhesive Base

KS 1069

Elastis Adhesive Dasar
Adhesivo Con Base Elastica
Elastic Base Adesiva
Elastic base adh sive

FEATURES

- Used for Fixation of I.V Cannula / Catheter
- Made from high-strength and extra breathable elastic fabric
- Uniformly coated with hypoallergenic, anti-fungal, non-irritant adhesive



➤ Sterile packed in peelable pouch pack.

STERILE

➤ I.V Transparent

KS 1070

I.V Transparan
I.V Transparente
i.v Transparente
i.v Transparent

FEATURES

- Used for Fixation of I.V Cannula / Catheter
- Made from Thin Polyurethane film and available in two variants
-Transparent and Non-woven
- Breathable and water proof, with reinforced notches that reduce mechanical stress from heavy catheters
- Thin film ensures easy fixation of devices with contoured surfaces



I.V. Infusion Set with Airvent

KS 1071

- I.V. Set Infusi dengan Airvent
- I.V. Set de infusión con airvent
- I.V. Conjunto de infusão com Airvent
- I.V. Set de perfusion avec Airvent

FEATURES

- IV Infusion Set with an airvent integrated above the drip chamber to allow adequate pressure inside the chamber for smooth, uninterrupted fluid flow.
- The airvent is fitted with a bacterial filter with snap cap to keep out the air borne pathogens.
- Transparent & Flexible drip chamber to facilitate rapid adjustment of fluid. Smooth roller clamp facilitates easy, safe control and adjustment of fluid rates
- Soft and Kink resistant PVC tubing ensure uniform flow rate, 150 cm length
- Sharp, sleek spike with vents to allow effortless penetration into the IV solution bottle.
- Moulded, self-sealing latex bulb at distal end for intermittent medication.
- Hypodermic needle attached to the male luer to prevent contamination.
- Soft, kink resistant tubing.
- Disposable, sterile and non-pyrogenic.
- Needle Size 21 G
- Approximate 20 drop / ml



► Sterile packed in peelable pouch pack.

STERILE

Measured Volume Burette Set

With "Y" Injection Port and Luer Lock

KS 1072

- Set Buret Volume Terukur
- Conjunto de bureta de volumen medido
- Conjunto de Buretas de Volume Medido
- Ensemble de burette à volume mesuré

FEATURES

- IV infusion set with a well-marked drip chamber for precise fluid administration.
- Flexible graduated burette with marking upto 100/150 ml.
- Sharp bevelled, smooth spike for easier penetration.
- Air inlet with bacterial filter and injection port.
- Disc filter inside the drip chamber for adequate filtration.
- C-clamp for regulating the rate of fluid infusion to the patient.
- Soft roller clamp for regulating the rate of fluid infusion to the patient.
- Kink resistant tubing 150cm with self-sealing latex bulb for extra medication.
- With Y injection port for intermittent medication and luer lock.
- Large drip chamber capacity 150 ml and small clamp.
- Needle Size 23G x 1.0"
- 60 Drops per mil Distilled Water.



► Sterile packed in peelable pouch pack.

STERILE

➤ Infusion Set with Filter

KS 1073

Set Infus dengan Filter
Set de infusión con filtro
Set de perfusion avec filtre
Conjunto de infusão com filtro

FEATURES

- IV Infusion Set with an airvent integrated above the drip chamber to allow adequate pressure inside the chamber for smooth, uninterrupted fluid flow.
- The airvent is fitted with a bacterial filter with snap cap to keep out the air borne pathogens.
- Secondary filter 0.2u to retain Bacteria / fungi and to remove entrapped air Disk filter to eliminate any particular matter & entrapped air bubble
- Transparent & Flexible drip chamber to facilitate rapid adjustment of fluid.
- Smooth roller clamp facilitates easy, safe control and adjustment of fluid rates
- Soft and Kink resistant PVC tubing ensure uniform flow rate, 150 cm length
- Sharp, sleek spike with vents to allow effortless penetration into the IV solution bottle.
- Moulded, self-sealing latex bulb at distal end for intermittent medication.
- Hypodermic needle attached to the male luer to prevent contamination.
- Soft, kink resistant tubing.
- Disposable, sterile and non-pyrogenic.
- Needle Size 21 G
- Approximate 20 drop / ml



STERILE

➤ Blood Transfusion Set Single / Double Chamber

KS 1074

Transfusi Darah Mengatur Kamar Tunggal / Ganda
Conjunto de transfusión de sangre de cámara simple / doble
Kit de transfusion sanguine simple / double
Conjunto de transfusão de sangue simples / câmara dupla

FEATURES

- Transfusion set specially designed to transfuse blood
- Strong, sleek and sharp spike with airvent to puncture blood bag outlets
- Drip chamber with filter mesh to keep out blood clots from entering the patient's circulation.
- Regulator clamp to control transfusion rate
- Kink resistant, soft, PVC tubing
- Injection site to administer medication
- Available in Single-drip/Double-drip chamber variants
- Options of Latex bulb site/Y-injection site Luer-slip/ Male-rotator variants
- Disposable, sterile and non-pyrogenic



3 Way Stop Cock With Extension Line

KS 1073	3 Way Stop Cock With Extension Line
	3 Llave De Paso Con La Línea De Extensión
	3 Way Parar Galo Com a extensão da linha
	3 Way robinet d'arrêt avec l'extension Ligne

FEATURES

- Used for pressure infusion and blood pressure monitoring
- Made from non-toxic, medical grade PVC
- Has two female luer lock ports and one male luer lock port, with rotatory lock for safe connection to other devices
- Three arm handle for smooth and easy operations
- Minimum residual fluid volume
- Optional- 25/50/100/200cm lengths available in extension line



› Sterile, individually packed in peelable pouch pack.

STERILE

3 Way Stop Cock

KS 10774	3 Way Berhenti Cock
	3 Way Llave De Paso
	3 Way Parar Galo
	3 Way robinet d'arrêt

FEATURES

- Used for pressure infusion and blood pressure monitoring
- Made from non-toxic, medical grade PVC
- Has two female luer lock ports and one male luer lock port, with rotatory lock for safe connection to other devices
- Three arm handle for smooth and easy operations
- Minimum residual fluid volume



› Sterile, individually packed in peelable pouch pack.

STERILE

P.E. Pressure Monitoring Line Spiral

Spiral

KS 1033

P.E. Pressure Monitoring Line Spiral
P.E. Línea Espiral De Monitoreo De Presión
P.E. Linha de Monitorização da Pressão
P.E. Surveillance de la pression de ligne

FEATURES

- Used for high pressure monitoring, and connection between syringe infusion pump and patient
- Made from Kink resistant PE Tubing
- Provided with universal male luer lock and female luer lock at either end to ensure securely fit all standard equipment
- Internal Diameter 1.0 mm, Outer Diameter 2.0 mm



› Sterile, Individually packed in peelable pouch & blister pack.

STERILE

Microporous Surgical Tape

KS 1075 S

- Tape mikroporous
- cinta microporosa
- Tape microporosa
- Tape microporeux

FEATURES

- Used for fixation of dressing, tubing, osteomy pouches etc in post-surgical care
- Made from Microporous fabric and specially coated with a broad spectrum Anti-Microbial agent for enhanced infection control
- A one-of-its-kind product, it minimizes the risk of the tape itself being a carrier of infection
- Hypoallergenic, non-toxic, non-irritant uniform adhesive coating
- Manufactured as per British Pharmacopia (BP) guidelines



SIZES

1/2", 1", 2" x 9.1mtr

3" x 9.1mtr

1/2", 1", 2" x 5.0mtr

3" x 5.0mtr

Surgical Tapes / Sports Tapes / Zinc Oxide Plaster


Uniformly Extra Porous

KS 1078

- Zinc Oxide Plester Seragam ekstra Porous
- Óxido de Zinc Yeso Uniformemente extra Poroso
- Óxido de Zinco Gesso Porous uniformemente extra
- Zinc Oxide Plaster Uniformément supplémentaire poreux

FEATURES

- Used for dressing and general medical applications
- Made from cotton blended fabric
- Uniformly coated with hypoallergenic, non-toxic, non-irritant adhesive
- Easy to peel off without causing discomfort or skin irritation
- Thermo sensitive adhesive activates better when it comes in contact with warm human skin
- Manufactured as per British Pharmacopia (BP) United State Pharmaeopia (USP) guidelines



SIZES

Width	Length
2.5cm	5mtr.
5.0cm	5mtr.
7.5cm	5mtr.
10.0cm	5mtr.

Transparent Perforated Plastic Film Tape

KS 1077 C

Transparan Tape berlubang Plastik Film
Cinta de película plástica transparente perforada
Fita de filme plástico transparente perfurada
Ruban de film plastique transparent perforé

FEATURES

- Made from perforated plastic film for good breathability
- Uniformly coated with hypoallergenic adhesive
- Strong adhesion water resistant
- Painless removal
- Perfect bidirectional tearing



SIZES

1/2", 1", 2" x 9.1mtr

3" x 9.1mtr

1/2", 1", 2" x 5.0mtr

3" x 5.0mtr

Sterile Absorbent Gauze Swab for Dressing

KS 1079 A

Swab Kasa Penyerap Steril untuk Berpakaian
Esponja de gasa absorbente estéril para vestirse
Cotonete absorvente estéril da gaze para vestir-se
Tampon de gaze absorbant stérile pour la vinaigrette

FEATURES

- Usage: Used for Dressing, General Wound Care, and Wound Cleaning.
- Size: 5cm x 5cm, 7.5cm x 7.5cm, 10cm x 10cm
- Material: 100% Cotton Swabs
- Strength: Available in 4ply, 8ply, 12ply And 16ply
- Features: High Absorbency Level. With / Without X-Ray Detectable Thread.
With / Without Edge Folded In Optical
- Sizes (cm) : 25 x 25, 30 x 30 & 35 x 35.



Available in Sterile/Non Sterile pack.

STERILE

► Surgical Pad / Gauze & Cotton Tissue

KS 1080

Surgical Pad/Gauze & Cotton Tissue
Pad Quirúrgica/Gasa Y Tejidos De Algodón
Pad cirúrgico / Gauze & Cotton Tissue
Pad chirurgicale / Gaze & Coton Tissu



FEATURES

- Used for dressing and general medical applications
- Made from fine absorbent cotton covered with a layer of non-woven gauze fabric
- Manufactured as per British Pharmacopia (BP)/USP guidelines
- Sizes (cm) : 6 x 6, 10 x 10 & 10 x 20.

► Elastic Adhesive Bandage BP

KS 1081

Elastis Adhesive Bandage BP
Elástico Esparadrapo BP
Elastic Curativo BP
Elastic Bandage BP



FEATURES

- Made from superior highly quality woven cloth with special porous adhesive permeable to air and moisture
- Used to support a weak joints as a stretch compression bandage in clinical and veterinary setting
- 100% Latex Free coating, No latex allergies and skin irritation
- Available in different size options
- Sizes: 6 cm x 4/6 mtr, 7.5 cm x 4/6 mtr, 8 cm x 4/6 mtr, 8 cm x 4/6 mtr, 10 cm x 4/6 mtr, 15 cm x 4/6 mtr

► Kay Crepe Bandages

KS 1076

Perban Kay Crepe
Vendas Kay Crepe
Kay Crepe Bandages
Bandages Kay Crêpe

FEATURES

- Reduces swelling
- Secures wound dressing
- Supports strains and sprains
- Kay Crepe Bandages are strong support crepes for medium to long-term bandaging requirements.
- They are ideal for securing wound dressing, supporting strains, and providing compression to alleviate swelling.
- Sizes: 6 cm x 4/6 mtr, 7.5 cm x 4/6 mtr, 8 cm x 4/6 mtr, 8 cm x 4/6 mtr, 10 cm x 4/6 mtr, 15 cm x 4/6 mtr 36



Manual Vacuum Aspiration (MVA Kit)

Double Valve / Single Valve

KS 1082

Aspirasi Vacuum Manual (MVA Kit)
La Aspiración Manual Endouterina (Ameu Kit)
Manual de Aspiração a vácuo (MVA Kit)
Aspiration manuelle (MVA Kit)

Double Valve MVA Syringe



Single Valve MVA Syringe

FEATURES

- Used for Manual Aspiration procedure (Medical termination of pregnancy)
- Made from medical grade polypropylene
- Kit consists of MTP Syringe, Karman Cannula, lubricating oil
- The MTP Syringe provides high vacuum suction during abortion procedure
- Smooth cannula with rounded tip and smooth eyes for atraumatic cannulation
- Available in single valve and double valve syringe
- Karman Cannula available separately as well in size 4mm - 12mm

➤ Sterile, individually packed in peelable pouch pack.

STERILE

Cytobrush / Endocervical Brush

KS 1083

Cytobrush / Endoserviks Sikat
Citocepillo/Cepillo Endocervical
Cytobrush / Endocervical Escova
Cytobrosse / Endocervical Brush

FEATURES

- Used for Cervical Papanicolaou (Pap) smear screening
- Consists of Brush and the Wooden spoon Strip



➤ Sterile, individually packed in peelable pouch pack.

STERILE

► Infant Mucus Extractor

KS 1084

Infant Mucus Extractor
Aspirador Infantil De Mucosidades
Infant Muco Extractor
Infant mucus Extractor

FEATURES

- Used for aspiration of secretion from oropharynx in newly born babies
- Made from non-toxic, medical grade PVC
- Transparent chamber for visual examination of the content
- Round open tip for atraumatic intubation, and proximal end with funnel connector for easy connection to the suction equipment
- Optional: available with filter and without filter models



SIZES

Chamber Capacity	25ml
Catheter Size	10 FG

► Sterile, individually packed in peelable pouch pack.

STERILE

► Umbilical Cord Clamp

KS 1085

Umbilical Cord Clamp
Abrazadera Cordón Umbilical
Cordão Umbilical Braçadeira
Pince cordon ombilic

FEATURES

- Used for clamping the umbilical cord of newly born babies
- Made from non-toxic, medical grade PVC
- Permanent locking system to prevent clamp reopening
- Finger grip to facilitate easy holding and locking of the clamp



SIZES

1/2", 1", 2" x 9.1mtr
3" x 9.1mtr
1/2", 1", 2" x 5.0mtr
3" x 5.0mtr

► Sterile, individually packed in peelable pouch pack.

STERILE

Umbilical Catheter

KS 1086


Kateter Umbilical
Catéter Umbilical
cateter umbilical
cathéter ombilical

FEATURES

- Used to access the umbilical vessels of pre-mature or newly - born babies
- Made from non-toxic, medical grade PVC
- Smooth round tip for atraumatic insertion and well-finished surface facilitate smooth passage in the vein
- Markings at every cm to ascertain depth of insertion, radio opaque line ascertain position during scans
- Proximal end provided with Luer female connector for easy connection to the I.V. therapy equipment

SIZES	
Total Length	35cm

4 White 5 Grey 6 Green 8 Blue 10 Black



➤ Sterile, individually packed in peelable pouch pack.

STERILE

Delivery Kits - All About Safe Child Birth

Normal and Caesarian

KS 1087

Kits Pengiriman Semua tentang pengaman pada saat kelahiran Anak
Kits De Entrega Todo sobre Nacimiento Niño Seguro
Kits de entrega - Tudo sobre Nascimento Criança Segura
Kits de livraison - All About Child Safe Naissance

FEATURES

- Complete solution for clinicians in normal and Caesarian delivery procedures
- Ensures safe child birth in case of unassisted deliveries in rural areas
- Standardized / Customized kits available

CONTENT OPTIONS FOR CUSTOMIZABLE KITS:

- Gauze Swab
- Abdominal Sponge (X-Ray op.)
- Umbilical Cord Clamp
- Umbilical Cord Ties
- Umbilical Tabular band (Elastic)
- Umbilical Scissor
- Sutures
- Forceps
- Needle Holder
- Mucus Extractor with Adopter
- Syringe (20 ml)
- Baby Wipe (60x60)
- Baby Diaper
- Baby Blanket
- Baby Garment
- Feeding Bottle
- Maternity Gown
- Maternity Pads
- Maternity Mesh Panty
- Surgeon Caps
- Surgeon Gown
- Surgeon Gloves
- Lignocaine Gel
- Povidone Iodine Solutions
- Plastic Tray
- Applicator Stick
- Perineal Sheet with Legging
- Maternity Drape Sheet
- Under Pad with Collecting Bag
- Instruction Sheet
- Measuring Tape
- Dimension Detail Sheet
- Pen
- ID Band (Pediatric)
- ID Band Adult
- Waste Bag



Colostomy Bag

KS 1016

Tas Kolostomi
Bolsa de colostomía
Sac de colostomie
Saco de colostomia

FEATURES

- Colostomy Bag is used for closed colostomy procedures.
- Manufactured from clinical grade PVC sheet.
- Coated with Hypo-allergic adhesive to prevent skin maceration.
- Hole diameter and bag can be resized as per customer's requirement.

SIZES	
Bag Size	Hole Diameter
250mm x 140 mm	20, 30, 40, 50mm



➤ Pack of 10

Colostomy Kit

KS 1017

Kit Kolostomi
Kit de colostomía
Kit de colostomie
Kit de Colostomia

FEATURES

- The Colostomy Kit gives perfect relief and comfort during prolonged use
- It is light weight, flexible and the belting arrangement fits the body well
- 'U' form padding is provided between the sleeve and main body for better comfort
- 100 poly spare bag provided.




➤ Individually Packed in a box

EUROVEIN Blood Collection Tubes Vacuum / Non Vacuum

KS 1091 EUROVEIN Tabung Pengumpul Darah Vacuum / Non Vacuum
 Tubos de extracción de sangre EUROVEIN Vacío / No vacío
 Tubes de prélèvement sanguin EUROVEIN avec ou sans vide
 Tubos de coleta de sangue EUROVEIN vácuo / sem vácuo

EUROVEIN Blood Collection tubes are manufactured under the strict surveillance of the experts using advanced machines and equipment. The use of quality plastic enables to preserve the stored blood without any changes to its size of the leukocyte and maintain invariable erythrocyte. EUROVEIN range of vacuum tubes are widely acclaimed for their features like light weight, excellent finishing and durability. EUROVEIN Blood Collection Tubes are available with Vacuum & Without Vacuum



WE OFFER THE FOLLOWING BLOOD COLLECTION TUBES

- EDTA Tube
- Heparin Tube
- Fluoride Tube
- No Additive Tube
- Coagulation Tube
- Clot Activator Tube
- Gel And Clot Activator Tube.

SIZES	
13 x 75 mm	4ml
13 x 100 mm	5ml
16 x 100 mm	10ml

Available in Sterile pack.

STERILE

Urine/Stool Specimen Container

KS 1088

Wadah Urin
 Recipiente de orina
 Recipiente de Urina
 Récipient d'urine

FEATURES

- Urine/Stool Specimen Container is used for the Collection of Urine /Stool Sample for the laboratory examination
- Made from medical grade Translucent Polypropylene or clear polystyrene

Urine Specimen Container

KS 1089

FEATURES

- Urine Specimen Container provided with easy to read graduation helps to determine the Volume.
- Provided with airtight Lid for safe handling of sample
- Available in 30ml and 70ml capacity.



➤ Sterile, individually packed in tear pouch pack.

STERILE

ECG Gel / Ultrasound Gel

KS 1090

ECG Gel / Gel Ultrasound
 Gel ECG / Gel Ultrasonido
 Gel ECG / Gel Ultra-som
 ECG Gel / Gel à Ultrasons

FEATURES

- Suitable for Short & Long Term Monitoring
- Non-irritating to Skin
- Neutral pH
- No Staining
- High Conductivity
- For external use only




Blood Lancet

KS 1063

Lancet darah
Lanceta de sangre
Lancette de sang
Lanceta de Sangue

FEATURES

- Blood Lancet used as a manual blade to make a limited skin incision as for collection of Blood.
- Made from medical grade Stainless Steel and polymer.
- Ultrasonic cleaned.
- Super sharp tip for the instant incision.



› Sterile, individually packed in peelable pouch pack.

STERILE

Alcohol Pad

KS 1065

Alkohol Pad
Almohadilla de alcohol
Tampon d'alcool
Almofada de álcool

FEATURES

- Alcohol pads are nonwoven strip saturated with 70% (v/v) Isopropyl Alcohol.



› Individually packed in paper aluminum foil laminate pouch with good barrier properties, the product has good self life.



Medical Disposables

Catalogue 2022

Representative:



Kay & Company

Head Office:

25, Netaji Subhash Marg, New Delhi - 110002. INDIA
Ph. : +91 11 43195600 - 20 Lines
Email : info@kay-co.com, sales@kaycoindia.com

Works:

B-316, Okhla Industrial Area, Phase-1, New Delhi - 110020. INDIA
Ph. / Fax : +91 11 41613432 / 33
Email : info@kay-co.com

Website

www.kaycoindia.com
www.medicaldisposabledevices.com